

# Comfort

Love. Care.  
Comfort.



Product Catalog

2022

# ABOUT COMPANY

We are a Polish company that for the last 35 years has been producing specialistic rehabilitation strollers Comfort. They are designed for children and adults suffering from cerebral palsy (CP), spiral muscular atrophy (SMA) or other afflictions causing palsies or limbs paresis

## ONE OF OUR BASIC PRINCIPLES IS TO MEET INDIVIDUAL EXPECTATIONS OF OUR CLIENTS.

While creating Comfort strollers we focus on their optimal influence on the patients. Bearing in mind how important the vestibular and proprioceptive systems are for the patient's growth, we pay particular attention to keeping the good body posture in a stroller.

## OUR STROLLERS ARE REGISTERED IN THE OFFICE FOR REGISTRATION OF MEDICINAL PRODUCTS, MEDICAL DEVICES AND BIOCIDAL PRODUCTS.

In 2005 we implemented the international system of standardisation ISO 1901:2000, which results in meeting the highest standards of production and customer service.

Our aim is to provide our clients with comfort in every sphere of their lives. That is why we actively support organizations caring for the disabled.

In compliance with the CSR (Corporate Social Responsibility) strategy introduced in our company, we get involved into numerous projects aiming at securing support and coming to help individuals with multiple and complex disabilities.





# FEATURES OF COMFORT STOLLERS



## Comfort

All reviews concerning Comfort wheelchairs show that they are supremely comfortable for children, teenagers as well as adults. Comfort wheelchairs are very easy to use.



## Reliability

We have 40 years of experience in producing specialistic wheelchairs. Hence, we can offer a 25-year guarantee for the wheelchair frame.

Comfort wheelchairs are reliability because we never use plastic elements.



## Support

Configurator of the Comfort wheelchair will help you to choose adequate equipment.

We can recomend solutions for the maximum subsidy. We can service your product both under/out of warranty within few days.



## Individuality

Indyvidual design of Comfort wheelchair may include shelf for respirator, profiled upholstery under the hump or special leg straps.



## Safety

The construction and balance of Comfort wheelchair guarantee stability and safety of its user.

Tube tyres and two-point brake used in Comfort wheelchair guarantee both a safe stop of the wheelchair.



## Practicality

Optimal size enables you to fold it and place in a small car. Swivel wheels make it possible to move both in small spaces and in the field. Recumbent and semirecumbent option available.

Double wheelchair seat upholstery enables using it even when the outer layer needs to be cleaned. Construction of comfort strollers enable the child to sit at the table. You can choose from over a dozen of useful elements.

Comfort wheelchair can be adjusted to be used together with a stairclimber TGR Jolly Junior. Comfort stroller ensure the body stability due to the optimal seat profile and foam pads, which have a positive influence on pripreception.

Comfort wheelchair can be used together with stabilizing pillows Stabilo.



# SPECIAL COMFORT WHEELCHAIR FOR CHILDREN

Comfortable, safe, reliable, functional, fully ergonomic, a special wheelchair for children ensuring comfort and mobility every day. Perfect as a rest and stroller to move around the city, in the field and indoors



## DEEP ADJUSTABLE HEADREST

It stabilizes the correct arrangement of the neck and head with muscle laxity of the nape or with permanent degenerative changes.

## SPECIAL SEAT PROFILE

It guarantees stabilization and very high comfort even during many hours spent in a stroller. It will act positive for proprioception.

## STABILIZING PADS

-

## ABDUCTION WEDGE

Drains the lower limbs, gently correcting them knee alignment without undue pressure point.

## ADJUSTABLE FOOTREST

Easy adjustment of the footrest position horizontally and vertically depending on the needs user.

## SWIVEL FRONT WHEELS

They make it easier to drive and maneuver the stroller, especially in narrow and confined spaces and in the field.

## BACKREST WITH THE OPTION LYING

Smooth adjustment of the backrest angle to the position lying as an option.

## COMFORTABLE BELTS

Belts to stabilize the lower limbs, pelvis, chest depending on user needs.

## UPHOLSTERY

Made of easy-to-clean material.

## CONVENIENT DRIVEWAY TO THE TABLE

The construction of the stroller allows the child to sit down conveniently at the table (for selected models).

## PNEUMATIC WHEELS

Solid wheels used in the smallest trolleys guarantee full maintenance-free.



### Available colors





# SPECIAL COMFORT WHEELCHAIR FOR ADULTS

Comfortable, safe, reliable, functional, fully ergonomic, special wheelchair for adults providing comfort and mobility every day. Perfect as a rest and stroller to move around around the city, in the field and indoors



## DEEP ADJUSTABLE HEADREST

It stabilizes the correct arrangement of the neck and head with muscle laxity of the nape or with permanent degenerative changes.

## SPECIAL SEAT PROFILE

It guarantees stabilization and very high comfort even during many hours spent in a stroller. It will act positive for proprioception.

## STABILIZING PADS

-

## ABDUCTION WEDGE

Drains the lower limbs, gently correcting them horizontally and vertically depending on the needs user.

## ADJUSTABLE FOOTREST

Easy adjustment of the footrest position horizontally and vertically depending on the needs user.

## SWIVEL FRONT WHEELS

They make it easier to drive and maneuver the stroller, especially in narrow and confined spaces and in the field.

## BACKREST WITH THE OPTION LYING

Smooth adjustment of the backrest angle to the position lying as an option.

## COMFORTABLE BELTS

Belts to stabilize the lower limbs, pelvis, chest depending on user needs.

## UPHOLSTERY

Made of easy-to-clean material.

## CONVENIENT DRIVEWAY TO THE TABLE

The construction of the stroller allows the child to sit down conveniently at the table (for selected models).

## PNEUMATIC WHEELS

Solid wheels used in the smallest trolleys guarantee full maintenance-free.



## Available colors




# USERS AND SIZES

	Wheelchair for children				Wheelchair for adults		
	Baby [1]	Size [3]	Size [4]	Size [5]	Size [6]	Size [6+]	Size [7]
<b>The purpose of the wheelchair</b>							
Age of the patient	0 - 3	3 - 6	6 - 12	6 - 14	from 14	from 14	from 14
Maximum patient weight [kg]	30	35	50	50	75	75	90
Patient height [cm]	110	140	145	155	150 - 165	160 - 180	160 - 180
<b>Wheelchair parameters</b>							
Seat width [cm]	33,5	35	38	38	41	41	44,5
Seat length/depth [cm]	27	33	37	41	45	48	48
Backrest height [cm]	55	66,5	70	76,5	78	83	83
Maximum opening angle against the backrest [°] standard position	up to 120°	up to 135°	up to 130°	up to 135°	up to 120°	up to 120°	up to 120°
Footrest height adjustment range [cm]	23 - 26	27,5 - 33,5	31,5 - 42	32 - 42,5	36 - 46	36 - 46	36,5 - 46,5
Max footrest opening angle in relation to the seat [°]	110°-155°	110° -145°	100° -135°	95° -135°	100° - 140°	95° - 150°	90° -130°
Approximate weight of trolley [kg]	15	21	23	23	26,2	26,5	28
Stroller dimensions after folding (length x width x height) [cm]	97,5 x 55 x 50,5	89,5 x 62 x 55	94 x 64,5 x 51,5	98 x 64,5 x 55	112,5 x 69,5 x 56	122 x 69,5 x 55	122 x 70 x 54
Stroller dimensions (length x width x height) [cm]	111 x 55 x 92	123 x 62 x 96	121,5 x 64,5 x 94,5	122,5 x 64,5 x 94,5	141 x 69,5 x 95	143 x 69,5 x 95	148 x 70 x 95





# EQUIPMENT



**Fixed wheels 12,5" x 2,25**  
(Cat. no.: ST)



**Seatbelts**  
(Cat. no.: ACPASBEZPR)



**Double upholstery**  
(Cat. no.: ACTAPW)



**Barrier**  
(Cat. no.: ACBAR)



**Inter-body wedge**  
(abductive)  
(Cat. no.: ACKLIN)




**Bag for the stroller**  
(Cat. no.: ACTORPWW)



**Headrest**  
(Cat. no.: ACZAG)


# FUNCTIONAL SOLUTIONS



**Full swivel wheels**  
(Cat. no.: SK1)



**Adjustable handle**  
(Cat. no.: RŁ)



**Backrest adapted to the hump**  
(Cat. no.: ACTAPG)



**Anti-tilt protection**  
(Cat. no.: ACUPP)



**Backrest tilted to a lying or reclining position [°]**  
(Cat. no.: OL)



**Raised armrests**  
(Cat. no.: PP)

# ADDITIONAL ACCESSORIES




**A booth with an apron**  
(Cat. no.: BUDA)



**Thigh belts**  
(Cat. no.: ACPASU)



**Cape**  
(Cat. no.: ACPEL)



**Side cushions**  
(Cat. no.: ACPODUB)



**Sleeping bag**  
(Cat. no.: ACŚPI)



**Lumbar supports - pelottes**  
(Cat. no.: ACPODPL)



**Umbrella**  
(Cat. no.: PARASOL.KOL.)



**Lumbar pillow**  
(Cat. no.: ACPODUL)



**Table**  
(Cat. no.: ACSTOL)




**Seat cushion**  
(Cat. no.: ACPODUS)



**Handle bag**  
(Cat. no.: TORBA.KOL.)




**Outer upholstery**  
(Cat. no.: ACTAPW)




**Hand protectors**  
(Cat. no.: ACOCHRAC)



**Footrest protector**  
(Cat. no.: ACOCHPOD)



**Calf belts**  
(Cat. no.: ACPASŁ)



**Ankle belts**  
(Cat. no.: ACPASK)



**Foot belts**  
(Cat. no.: ACPASS)



# ADDITIONAL INFORMATION

We invite you to use the electronic media library and to follow our social media.



## YouTube

Channel presenting Comfort film materials.  
[youtube.com/channel/UCj8kE1tWS87ry4be080OpHQ](https://youtube.com/channel/UCj8kE1tWS87ry4be080OpHQ)



## Facebook

Brand profile informing about new products, current activities of the company and charity activities.  
[facebook.com/Comfort-w%C3%B3zki-rehabilitacyjne-105814814917073](https://facebook.com/Comfort-w%C3%B3zki-rehabilitacyjne-105814814917073)



## Instagram

An image profile showing the use of prams, the use of accessories and functional solutions.  
[instagram.com/comfortpolska](https://instagram.com/comfortpolska)




## Digital Content

Electronic media library where you can find and download for your digital presentations: photos, videos, user descriptions and opinions.  
[comfort-reha.pl/content/](https://comfort-reha.pl/content/)




# MEASUREMENT TABLE

How to correctly measure the User? In order to choose the right one the size of the trolley, the measurement of the User should be made.



The most important are 5 dimensions:

- 1 The length of the body from the seat to the top of the head
- 2 The length of the body from the seat to the shoulders
- 3 The length of the body from the back to the knee well
- 4 Shin length (from the knee to the sole)
- 5 Hips (distance between the widest places)



If there is a hump, it is necessary accurate dimensioning.

- 1 Length of the lumbar spine
- 2 The length of the hump from the shoulders to the end of its occurrence
- 3 The length of the body from the seat to the top of the head

Dodatkowo należy wskazać czy garb jest lewostronny czy prawostronny.





Love. Care. Comfort.

**Comfort**

**Factory of Rehabilitation Equipment**

ul. Krzywa 23

60-118 Poznań, Poland

phone +48 61 863 85 61

biuro@comfort-reha.pl

**[www.comfort-reha.pl](http://www.comfort-reha.pl)**