

Comfort[®]

PRODUCT CATALOGUE

**Specialised respite
strollers for children
and adults**



Love. Care. Comfort.

This is a medical device. Use it in accordance with the instructions for use.

Care. Comfort. In practice.

Care for the user starts with a well-thought-out design. That is why Comfort strollers combine comfort, safety, functionality and an individual approach with more than 40 years of experience.



Reliability

More than 40 years of experience, together with our own production and design department in Poland, give us full control over the quality of every stroller.



Comfort

Supportive seating and adjustable headrests help provide comfort and stable support for the head and neck.



Support

We help select the right model and equipment, support our partners and provide efficient servicing.



Individual approach

We tailor the stroller not only to the user's height and weight, but also to their individual needs.



Safety

A foot brake, safety bar and abduction wedge help keep the user secure and reduce the risk of sliding out of the seat.



Functionality

Adjustable backrest and footrest settings make it easier to tailor the stroller to the user's individual needs.

We start with a conversation

At Comfort, we do not start with a ready-made solution. We start with a conversation - with the user, caregiver or our business partner. We want to understand everyday needs and challenges so that we can propose a solution that genuinely addresses them.



For more than 40 years, we have been developing our products with input from users, healthcare professionals and distributors.

Their feedback helps us improve our strollers and develop solutions shaped by everyday experience.

We combine our experience as a manufacturer with the local knowledge of our distributors, allowing us to stay close to customers across different markets. We also support NGOs and projects dedicated to helping people with disabilities and their caregivers.

Comfort

TRIP

Frequent travel often involves repeated transfers between the stroller and the vehicle seat. This takes time and can be physically demanding for both the user and caregiver.

Perhaps you already have a vehicle adapted for transporting people with disabilities, but the stroller you use has not been designed and tested for transporting the user while they remain in it?

Comfort TRIP brings two functions together in one solution: a specialised stroller for everyday use that is also designed for transporting the user in a vehicle adapted for people with disabilities.



Compliant with ISO 7176-19:2022



Dynamic crash testing confirms the stroller meets the performance requirements of ISO 7176-19:2022 under simulated frontal-impact conditions.

Four-point securement



Two securement points at the front and two at the rear allow the stroller to be securely restrained using an ISO 10542-1-compliant four-point tiedown system.

Versatile use

Suitable for vehicle transport as well as everyday use at home, at school and during rehabilitation.

Reinforced frame



Designed to withstand the forces acting on the stroller during braking and other road situations.

Key Features

Swivel front wheels make manoeuvring easier, the folding push handle reduces the stroller's overall length, making it easier to fit into a car, while the reclining backrest allows the user to lie down.

Technical specifications

Size	3	4	5	6	6+	7
USER						
Approximate user age [years]	3–6	6–12	6–14	od 14	od 14	od 14
Max. user weight [kg]	35	50	50	75	75	90
Max. user height [cm]	140	145	155	155–165	160–180	160–180
SEAT						
Seat width [cm]	35	38	38	41	41	44,5
Seat length [cm]	33	37	41	45	48	48
Backrest height [cm]	66,5	70	76,5	78	83	83
Max. backrest-to-seat angle [deg]	up to 150	up to 145	up to 150	up to 140	up to 140	up to 140
ADJUSTMENT						
Footrest height adjustment range [cm]	27,5–33,5	31,5–42	32–42,5	36–46	36–46	36,5–46,5
Footrest-to-seat angle [deg]	110–145	100–135	95–135	100–140	95–150	90–130
PARAMETERS						
Approximate stroller weight [kg]	24	26	26	29,2	29,5	31
Stroller dimensions when folded L x W x H [cm]	89,5 x 62 x 55	94 x 64,5 x 51,5	98 x 64,5 x 55	112,5 x 69,5 x 56	122 x 69,5 x 55	122 x 70 x 54
Stroller dimensions when unfolded L x W x H [cm]	123 x 62 x 96	121,5 x 64,5 x 94,5	122,5 x 64,5 x 94,5	141 x 69,5 x 95	143 x 69,5 x 95	148 x 70 x 95

Available colour options



Scan the QR code
to find out more



Comfort

Easy Vero

A few years ago, transferring the user from the stroller to a bed, sofa or bath may not have been particularly difficult. But as they have grown and gained weight, have transfers become more physically demanding?

Or does everyday care require easy, unrestricted access to the user that fixed side sections make difficult?

Comfort Easy Vero is a specialised rehabilitation stroller that combines the functionality, quality and durability associated with Comfort strollers with easily removable sides, making transfers and everyday care easier.



Open design



Removable sides provide easier access to the user, making transfers, repositioning and everyday care simpler.

Compatible with patient lifts



The removable sides provide easier access to the user and facilitate transfers with most types of patient lifts.

Quick-release abduction wedge



An abduction wedge with a quick-release mechanism for easy installation and removal.

Easy folding and transport

The stroller folds easily and fits into the boot of a smaller car.

Technical specifications

Size	3	4	5	6	6+	7	8
USER							
Approximate user age [years]	3–6	6–12	6–14	od 14	od 14	od 14	od 14
Max. user weight [kg]	35	50	50	75	75	90	130
Max. user height [cm]	140	145	155	150–165	160–180	160–180	160–180
SEAT							
Seat width [cm]	35	38	38	41	41	44,5	52,5
Seat length [cm]	33	37	41	45	48	48	48
Backrest height [cm]	66,5	70	76,5	78	83	83	83
Max. backrest-to-seat angle [deg]	up to 135–150	up to 130–145	up to 135–150	up to 120–140	up to 120–140	up to 120–140	up to 120–140
ADJUSTMENT							
Footrest height adjustment range [cm]	27,5–33,5	31,5–42	32–42,5	36–46	36–46	36,5–46,5	36,5–46,5
Footrest-to-seat angle [deg]	110–145	100–135	95–135	100–140	95–150	90–130	90–130
PARAMETERS							
Approximate stroller weight [kg]	22	24	24	27,2	27,5	29	33
Stroller dimensions when folded L x W x H [cm]	89,5 x 62 x 55	94 x 64,5 x 51,5	98 x 64,5 x 55	112,5 x 69,5 x 56	122 x 69,5 x 55	122 x 70 x 54	122 x 78 x 95
Stroller dimensions when unfolded L x W x H [cm]	123 x 62 x 96	121,5 x 64,5 x 94,5	122,5 x 64,5 x 94,5	141 x 69,5 x 95	143 x 69,5 x 95	148 x 70 x 95	148 x 78 x 95

Available colour options



Z04 Z16 Z18

Scan the QR code
to find out more



Specialistic rehabilitation stroller Comfort

Baby

Choosing your child's first specialised stroller is an important decision. As well as providing the right support, the caregiver's comfort also matters.

The stroller should be easy to manoeuvre and fold, fit comfortably into a car boot and be easy to keep clean. Are you looking for a lightweight, compact solution designed for the youngest users that also makes everyday care easier?

Comfort Baby is the smallest model in the Comfort range, designed for the youngest users who require appropriate support for the head, trunk and limbs. Its design helps provide the child with a stable and comfortable position, while a wide range of accessories allows the stroller to be tailored to the user's individual needs.



Compact design



The stroller folds easily and fits into the boot of a smaller car.

Adjustable push handle



Allows the push handle to be adjusted to suit the caregiver, making the stroller more comfortable to manoeuvre.

Backrest reclining to a lying position



Solid swivel wheels

Make manoeuvring easier indoors and during outings, and thanks to their solid tyres, eliminate the risk of punctures.

Easy-clean upholstery

The upholstery fabric makes everyday dirt easier to remove.

Allows the backrest to be reclined, providing the child with a comfortable position for rest.

Technical specifications

Size	1
USER	
Approximate user age [years]	0-3
Max. user weight [kg]	30
Max. user height [cm]	110
SEAT	
Seat width [cm]	33,5
Seat length [cm]	27
Backrest height [cm]	55
Max. backrest-to-seat angle [deg]	up to 130-150
ADJUSTMENT	
Footrest height adjustment range [cm]	23-26
Footrest-to-seat angle [deg]	110-155
PARAMETERS	
Approximate stroller weight [kg]	15
Stroller dimensions when folded L x W x H [cm]	97,5 x 55 x 50,5
Stroller dimensions when unfolded L x W x H [cm]	111 x 55 x 92

Available colour options



Scan the QR code
to find out more



Specialistic rehabilitation stroller Comfort

Vario

Everyday life takes users through different environments — at home, out and about, and on more demanding terrain. A stroller should provide a comfortable, stable position throughout, while also accommodating individual differences in body proportions and the level of postural support required.

Are you looking for a versatile stroller that can be tailored to the user's individual needs?

Comfort Vario is our most versatile stroller line, known for its functionality, build quality and durability.

Available in a wide range of sizes, it suits both smaller and larger users. Extensive personalisation options allow it to be tailored to the user's needs.



Contoured upholstery



Supports proprioception and stable positioning, helping to maintain comfort during prolonged use.

Double-layer upholstery



The removable top layer can be taken off for cleaning, while the fixed lower layer allows continued use of the stroller.

Wide range of personalisation options



A broad choice of accessories and custom-made solutions allows the stroller to be adapted to the needs of both the user and caregiver.

Position adjustment

Adjustable backrest and footrest settings make it easier to tailor the stroller to the user's individual needs.

Wide range of sizes

Available in a variety of sizes to accommodate both smaller and larger users.

Technical specifications

Size	3	4	5	6	6+	7
USER						
Approximate user age [years]	3–6	6–12	6–14	od 14	od 14	od 14
Max. user weight [kg]	35	50	50	75	75	90
Max. user height [cm]	140	145	155	150–165	160–180	160–180
SEAT						
Seat width [cm]	35	38	38	41	41	44,5
Seat length [cm]	33	37	41	45	48	48
Backrest height [cm]	66,5	70	76,5	78	83	83
Max. backrest-to-seat angle [deg]	up to 135–150	up to 130–145	up to 135–150	up to 120–140	up to 120–140	up to 120–140
ADJUSTMENT						
Footrest height adjustment range [cm]	27,5–33,5	31,5–42	32–42,5	36–46	36–46	36,5–46,5
Footrest-to-seat angle [deg]	110–145	100–135	95–135	100–140	95–150	90–130
PARAMETERS						
Approximate stroller weight [kg]	21	23	23	26,2	26,5	28
Stroller dimensions when folded L x W x H [cm]	89,5 x 62 x 55	94 x 64,5 x 51,5	98 x 64,5 x 55	112,5 x 69,5 x 56	122 x 69,5 x 55	122 x 70 x 54
Stroller dimensions when unfolded L x W x H [cm]	123 x 62 x 96	121,5 x 64,5 x 94,5	122,5 x 64,5 x 94,5	141 x 69,5 x 95	143 x 69,5 x 95	148 x 70 x 95

Available colour options



Scan the QR code
to find out more



Specialistic rehabilitation stroller Comfort

Spazio

A larger user needs adequate space and load capacity without having to compromise on stable support. Finding a passive stroller that combines all these features can often be a major challenge.

Are you looking for a specialised stroller that gives a larger user the space, comfort and stable positioning they need?

Comfort Spazio is the largest model in the Comfort family of specialised strollers. It is not simply an enlarged version of the Vario line — its design has been adapted for higher loads and the needs of larger users, while retaining effective postural support, practical functionality and extensive personalisation options.



Reinforced frame



Specially adapted for users weighing up to 130 kg.

Spacious seat



A seat width of 52.5 cm provides more space for users who need a larger seating area.

Contoured upholstery



Supports proprioception and stable positioning, helping to maintain comfort during prolonged use.

Body stabilisation

A deep, adjustable headrest, belts and additional support elements help keep the head, trunk and limbs in the correct position.

Wide range of personalisation options

A broad choice of accessories and custom-made solutions allows the stroller to be adapted to the needs of both the user and caregiver.

Technical specifications

Size	8
USER	
Approximate user age [years]	od 14
Max. user weight [kg]	130
Max. user height [cm]	160–180
SEAT	
Seat width [cm]	52,5
Seat length [cm]	48
Backrest height [cm]	83
Max. backrest-to-seat angle [deg]	up to 120–140
ADJUSTMENT	
Footrest height adjustment range [cm]	36,5–46,5
Footrest-to-seat angle [deg]	90–130
PARAMETERS	
Approximate stroller weight [kg]	32
Stroller dimensions when folded L x W x H [cm]	122 x 78 x 95
Stroller dimensions when unfolded L x W x H [cm]	148 x 78 x 95

Available colour options



Scan the QR code
to find out more



Comfort

Coala

People with cerebral palsy and other neurological conditions may experience severe spasticity, marked muscle flaccidity or profound hypotonia combined with a lack of head control.

Some users may also experience persistent en bloc body stiffness, particularly in the pelvis, lower-limb joints and spine. When these symptoms are persistent, assuming and maintaining a seated or semi-seated position may be impossible.

Are you looking for a solution that allows such a user to spend time away from bed while remaining in a horizontal position?

Comfort Coala is a specialised rehabilitation stroller designed for users who remain primarily in a horizontal position.



Adjustable backrest



The backrest can be adjusted within a range of 150°–180°, allowing the user's upper body to be raised and positioned to suit their needs.

Shelves for medical equipment



Space for essential equipment such as a ventilator or suction unit.

Removable seat



The seat can be removed from the frame, with integrated carrying handles making it easier to carry.

Key Features

Swivel front wheels make manoeuvring easier, while the folding push handle reduces the stroller's overall length for transport.

Stabilising legs

Fold-out legs add extra stability when the stroller is stationary.

Technical specifications

Size	9
USER	
Approximate user age [years]	-
Max. user weight [kg]	75
Max. user height [cm]	180
SEAT	
Seat width [cm]	45
Seat length [cm]	100
Backrest height [cm]	40
Max. backrest-to-seat angle [deg]	up to 180
ADJUSTMENT	
Footrest height adjustment range [cm]	-
Footrest-to-seat angle [deg]	-
PARAMETERS	
Approximate stroller weight [kg]	39
Stroller dimensions when folded L x W x H [cm]	150 x 58 x 20
Stroller dimensions when unfolded L x W x H [cm]	170 x 78 x 90

Available colour options



Z16



Scan the QR code
to find out more



Available colour options

Choose from a selection of upholstery colour options to suit individual preferences. Available colours vary by model – see the relevant product page for details.

Navy & blue [Z-01]



Navy [Z-04]



Navy with toucan pattern [Z-06]



Grey & turquoise [Z-09]



Grey & orange [Z-11]



Navy & red [Z-12]



Grey & red [Z-13]



Grey [Z-16]



Black [Z-18]



Accessories

Backrest adapted to accommodate a spinal prominence
[ACTAPG]



Suction unit shelf
[ACPPSS]



Reclining backrest for lying and semi-reclined positions
[OL]



Ventilator shelf
[ACPPRES]



Raised armrests
[PP]



Padded footrest protector
[ACPOCHOD]



Padded push handle protector
[ACOCHRAAC]



Calf straps
[ACPASL]



Foot straps
[ACPASS]



Thigh straps
[ACPASU]



Ankle straps
[ACPASK]



Lumbar supports – pads
[ACPELOT]



Lumbar cushion
[ACPODUL]



Seat cushion
[ACPODUS]



Side cushions
[ACPODUB]



Lateral supports
[UCHW.PEL]



Stroller tray
[ACSTOL]




Quick-release abduction wedge
[KLIN.AB.SZ.MON]



Accessories

Swivel front wheels
[SK1]



Adjustable push handle [RL]



Push-handle bag
[TORBA.KOL]




Helios 150 patient lift
[PHEL]




Anti-tip device
[ACUPP]



Sun umbrella
[PARASOL.KOL]



Stroller footmuff
[ACSPI]



Canopy with apron
[BUDA]



Rain cape
[ACPEL]




ComfortMemo orthopaedic cushion (size S, M, L)
[ACPODZPAM]




L-shaped stabilising cushion (size S, M, L)
[POD-L-STAB]




Pneumatic wheel with inner tube [ACKST4]




Tubeless wheel
[ACKST5]




Inner tube - size selectable [ACDET]




Tyre
[ACOP12CZ]




Wheel cap
[ACKAP5]




Tablet holder
[ACUCHTA]




Phone holder
[ACUCHTE]




IV pole holder
[UCH.KROPL]




Bottle holder
[ACUCHBI]




Registration plate
[AC.TAB.REJ]




Sunrise R20 electric add-on drive (size S, L) [DOSTAW.ELER20S/L]




Additional safety bar
[ACBAR]




Additional abduction wedge [ACKLIN]



Additional safety bar protector [ACOCBAR]



Additional four-point safety belts [ACPASBEZPR]



Additional top upholstery
[ACTAPW]



Additional headrest
[ACZAG]



Additional bag under the stroller [ACTORPWW]




Choose the right stroller size

To select the correct stroller size, measure the user according to the instructions below.




Measurements should be taken with the user seated in their usual position. Accurate measurements help ensure the correct stroller size and configuration.

Four measurements are required to select the appropriate stroller size:

1. Seated height – from the seat surface to the top of the head
2. Seat depth – from the back of the pelvis to the back of the knee
3. Lower leg length – from the back of the knee to the sole of the foot
4. Hip width



If the user has a spinal prominence, take the following additional measurements:



1. Vertical position of the apex
2. Horizontal position of the apex
3. Highest point of the prominence base
4. Length of the prominence
5. Width of the prominence
6. Vertical distance between the outermost points of the prominence
7. Height of the apex

WSR Comfort sp. z o.o.

ul. Głazowa 43
60-116 Poznań
tel. +48 61 863 85 61
biuro@comfort-reha.pl
www.comfort-reha.pl



Comfort special stroller. Intended for the mobility of children and adults with cerebral palsy and other conditions that prevent independent mobility. The stroller must be pushed by a caregiver.

Comfort Baby, Vario and Spazio do not constitute trade names within the meaning of Regulation (EU) 2017/745 on medical devices. Comfort Baby, Vario and Spazio are product brands used solely for marketing purposes.